

# MDHearingAid™

## The Doctor's Choice for Affordable Hearing Aids™




“As an Ear, Nose, and Throat physician, I see people on a daily basis that would benefit from hearing aids, but cannot afford the hundreds or thousands of dollars often needed for them. **MDHearingAid™** was designed to be the **first line of medical-grade hearing aids under \$200.**”

Dr. Sreek Cherukuri, MD  
Board-Certified ENT Physician  
Founder of MDHearingAid Inc.

## MDHearingAid Pro™

- **Our Flagship Model**— A Professional Behind-the-Ear hearing aid that is Audiologist-tested and Doctor-recommended
- #1 Rated hearing aid on Amazon.com
- 45 dB of gain. Best choice for Mild, Moderate, or Moderately-Severe hearing loss
- Fits *discreetly* behind the ear. Comes with various sizes of sized clear ear domes/tips to fit most adult ears
- Volume dial, 2-Channel tone control switch to decrease background noise and customize to predominantly flat or high-frequency hearing loss
- Provides amplification of the critical frequencies of the human voice (1000Hz - 4000Hz) without overly amplifying



**MDHearingAid PRO**

Visit us online at [www.MDHearingAid.com](http://www.MDHearingAid.com)  
or call our  
US-Based Call Center **1-888-670-HEAR** (4327)

**Always Open**  
24 Hours a day, Every  
Day, All Year

# MDHearingAid Max™



- The MDHearingAid MAX is a larger-sized basic hearing aid. It is a reasonable option for those requiring the service of a caregiver, and those with limited use of their hands.
- Receiver-in-the-ear (RITE) hearing aid design is easier to see and handle (visible receiver button sits in your ear).
- 2-Channel tone control
- 40 dB of gain for mild to moderate hearing loss
- Different sizes of ear domes to fit most adults
- General amplification hearing aid that will help hearing the television and family members.



**MDHearingAid MAX**

## Mail Order Form

Complete this form, and mail\* to:  
**MDHearingAid**, 917 W. Washington Blvd.,  
Suite 202, Chicago, IL 60607

Please select (Price includes shipping & handling):

- 1 MDHearingAid Pro \$185.94
- 2 MDHearingAid Pro \$329.99 **Save \$41.89**
- 1 MDHearingAid Max \$125.94
- 2 MDHearingAid Max \$205.94 **Save \$39.99**

International orders add \$25 for shipping/handling

Name: \_\_\_\_\_

Address: \_\_\_\_\_

City/State/Zip: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

PAYMENT METHOD:  Money Order

Visa  Mastercard  Discover  Amex

Card# \_\_\_\_\_

Expiration Date \_\_\_\_\_ CVV<sup>†</sup> \_\_\_\_\_

Signature (Required): \_\_\_\_\_

By signing above, I agree that I am over age 18 and I have been advised by MDHearingAid that the Food and Drug Administration and the State of Illinois have determined that my best interest would be served if I had a medical evaluation by a licensed physician, preferably a physician who specialized in diseases of the ear, before purchasing a hearing instrument; or a test by a licensed audiologist or licensed hearing instrument dispenser utilizing established procedures and instrumentation in the fitting of hearing instruments. I do not wish either a medical evaluation or test before purchasing a hearing instrument.

\* Delivery confirmation recommended for security

† 3 or 4 digit number after credit card number on back of most credit cards

**45 Day Risk-FREE Trial! - Money Back Guarantee**

Visit us online at [www.MDHearingAid.com](http://www.MDHearingAid.com)  
or call our  
US-Based Call Center **1-888-670-HEAR** (4327)

**Always Open**  
24 Hours a day, Every  
Day, All Year